Argon medicinal W21.8

Argon médical W21.8

Argon medico W21.8

Medical argon W21.8

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Versione 04

ENGLISH

Instructions for use – Information for the user

Read the entire instructions for use carefully because they contain important information. The instructions for use are an integral part of the medical device and must always be available to the user, e.g., the physician. Before working with the medical device, you must know and understand its function and handling.

The medical device Medical Argon consists of compressed argon gas. The compressed gas is filled in dedicated aluminium cylinders as the primary packaging, using a valve as the closure system.

Please note that a number of the points listed in this Instruction for Use may also pose a safety risk to the physician / user of Argon. Therefore, it is particularly important that the use of Argon must always be carried out by staff who have been trained with regard to the hazards associated with Argon.

For detailed information, see also chapters 5, 6, and 7.

Before any intervention, make sure that the product used is Medical Argon (marked CE 0123) by consulting the label.

Use a specific connection of type W 21,80 × 1/14 (Type 6 according to DIN 477) it is dangerous and forbidden to use similar or non-matching parts.

This instruction for use include:

1. Intended Purpose	65
2. Intended Clinical Benefits	65
3. Intended Use Environment	65
4. Indication, Dosage, Method and Duration of Use	66
5. Contraindications, Warnings, Precautions, Side-Effects	68
6. Precautions and Warnings for Transport and Storage	70
7. Safe Withdrawal	72
8. Cleaning and Disinfection	76
9. Maintenance, Servicing and Disposal	78
10. Further Information	79
11. Labelling	81

1. INTENDED PURPOSE

Medical argon gas is intended to be used as a medium for APC coagulators, whereby by connecting the gas in its primary packaging to the APC coagulator, the resulting argon plasma enables coagulation and the ablation of tissue anomalies during surgical procedures.

2. INTENDED CLINICAL BENEFITS

The intended effect of Medical Argon depends on its use in the procedure Argon Plasma Coagulation.

The decisive advantages of argon plasma coagulation are haemostasis and tissue revitalisation.

3. INTENDED USE ENVIRONMENT

Health care offices or institutions

No specific ambient conditions apply, other than those listed under section 6 "Precautions for Transport / Storage".

4. INDICATION, DOSAGE, METHOD AND DURATION OF USE

4.1 Medical Indication

Medical Argon is a medical device intended to supply the argon plasma coagulators, which are used during surgical procedures for haemostasis/coagulation or to ablate tissue in open surgery, endoscopy or laparoscopic surgery.

4.2 Intended Patient Population

No or limited use in patients with a clinical predisposition as described in the “Contraindication” and “Precaution” sections. The intended patient population is anyone who requires APC. Any limitations are specific to the APC procedure required. An individual risk assessment must be done by the medical doctor on a case-by-case basis.

4.3 Dosage and Duration of Use

The amount, speed and duration of application depend on the indication and are to be determined individually by the treating physicians, whenever possible, according to clinical guidelines.

Keep the flow rate to the lowest possible in case of danger of argon blowing into body cavities or gas embolism.

The duration of application is limited according to the classification rule «short-term» (i.e., up to a maximum of 24 hours).

4.4 Intended Users

Medical personnel (e.g., physicians)

Medical intervention may only be performed by physicians who have experience in the application of argon plasma coagulation.

Personnel for logistics of the cylinders

(e.g., technical personnel, nurses)

No demographic or physical restrictions when working according to the Linde Gas Schweiz safety data sheet.

4.5 Method and Application

Application should only be performed by physicians who have experience in the application of argon plasma coagulation. It should be performed with a state-of-the-art medical system. The user should always consult the manual of the APC device for procedure related risks.

5. CONTRAINDICATIONS, WARNINGS, PRECAUTIONS, SIDE-EFFECTS

- Should any incidents or side-effects occur during the use of Argon, which are not described in chapter 5.3, these must be reported to the manufacturer Linde Gas Schweiz.
- If you notice side-effects that are already described in chapter 5.3, you can also report them to the manufacturer Linde Gas Schweiz.
- If a serious incident occurs in relation to Argon, this serious incident must be reported immediately to the manufacturer Linde Gas Schweiz and to the competent authority of the Member State in which the user and/or the patient is established.

5.1 Contraindications

As Argon gas has to be used in combination with an APC device, there are no contraindications for the gas itself. For contraindications related to the APC device, the user should consult the device' information.

5.2 Warnings and Precautions

Comply with the instructions for use provided by the manufacturer of the device associated with the cylinder.

Danger of airway fire and intestinal explosion with other gases used in surgery or anaesthesia (e.g., oxygen, nitrous oxide). The argon plasma can be a source for thermal induction or ignition of those gases. The addition of oxygen and the use of nitrous oxide as an anesthetic must be strictly avoided.

5.3 Side-effects

Possible risks of complication associated with the use of gas with argon plasma coagulators:

- gas embolism
- pneumoperitoneum
- submucosal or subcutaneous emphysema
- intestinal gas distension

5.4 Interactions

Argon itself is an inert gas. For further information on argon plasma see section 5.2 “Warnings and precautions”.

5.5 Pregnancy and Lactation

The intended patient population is anyone who requires APC. Any limitations are specific to the APC procedure required. An individual risk assessment must be done by the medical doctor on a case-by-case basis.

5.6 Trafficability and the Operation of Machines

No special precautions are necessary.

6. PRECAUTIONS AND WARNINGS FOR TRANSPORT AND STORAGE

- The Argon is a compressed gas (200 bar at 15 °C) and it is supplied in pressurised cylinders.
 - Only store and transport gas cylinders with closed valves and with the protective device provided (e.g., protective cap).
 - Protect the valve from mechanical overstressing.
 - During storage, transportation, and use, restrain cylinders using an appropriate means (chains, hooks, etc.) to keep them in a vertical position and prevent them from falling.
 - Store at ambient temperature below 50 °C and do not expose to heat, may burst if heated.
 - Store in a clean, well-ventilated room. Gas is heavier than air. It can accumulate in closed rooms, especially on the floor or in low-lying areas.
 - Therefore, do not place cylinders in stairwells, corridors, passageways and recreation or consumption rooms.
 - Empty and full cylinders must be labelled and stored separately from each other.
- When storing empty cylinders ensure that the valves are kept shut and protected.

- The expiry date is printed on a separate batch label on the cylinder. Do not use the medical device after this date! The expiry date also applies to opened cylinders without restriction. The retest date of the regular ADR check of the receptacle is indicated on the neck of the cylinder; the retest date refers exclusively to the permission of transport by road and does not refer to the expiry date of the product.

7. SAFE WITHDRAWAL

7.1 General Safety Instructions

- Unless specifically stated otherwise, use cylinders with the valves at the top.
- Argon is heavier than air. At high concentrations there is a danger of asphyxiation. Use only in rooms with sufficient ventilation. Apply only by qualified personnel.

- An atmosphere in a closed room of less than 19.5 % oxygen can be dangerous; it can lead to loss of consciousness and later to asphyxia in confined space / small rooms ($< 15 \text{ m}^3$).
- In the event of argon gas escaping and it being impossible to evacuate it, the room must be left immediately. The room must not be entered until it is certain that the oxygen concentration has risen above 19.5 % again.
- In the event of gas escapes, close the leaking valve, ventilate the room well and evacuate it. Never use leaking cylinders and check that emergency measures have been taken.
- Check the cylinder before use to ensure that it is in perfect condition and suitable for the intended use.
- Remove the tamper-evident seal from the cylinder; ensure that the cylinder's fitting is clean and free of foreign bodies. Open the cylinder valve slowly until a

hissing sound can be heard, then close it immediately (this procedure purges any foreign bodies). Safety glasses are recommended.

- In case of any abnormalities do not use the product and contact the manufacturer.
- The cylinders are equipped with a residual pressure valve; the residual pressure is between 1.5 and 3.5 bar; below this pressure the cylinder no longer supplies gas.
- Do not lift cylinders by their valves.
- Never squeeze a cylinder into a holder, into which it is difficult to fit.
- Misuse of the cylinder, filling by the user or by third parties and decanting into other cylinders are not permitted. Do not transfer compressed gas from one cylinder to another.

7.2 Operating Argon Cylinders with a Pressure Reducer

- Before start of application make sure there is a sufficient residual amount to guarantee continuity of supply throughout the product delivery process. To do this, check the level by reading the pressure on the pressure gauge. An additional backup cylinder is recommended.
- Argon is not sterile and not free of particles. As this gas is not sterile, it must be filtered if necessary.

- Check the output connection of the cylinders for cleanliness before connecting the pressure reducer; clean any dirty connections with an oil- and grease-free cloth. Keep the interface between the cylinder and the pressure reducer clean. Check the condition of the seals.
- Only open cylinder with a pressure reducer connected. Before opening the cylinder valve, check the connected fittings (e.g., pressure reducer, flowmeter) for their closed condition. The pressure reducer must be relieved.
Open the cylinder valve slowly anticlockwise, never with force and as far as it will go. Tighten the pressure reducer/flowmeter preferably by hand to avoid damaging the seals.
- Never stand in front of the valve, but always position yourself on the side opposite the pressure reducer, behind the cylinder and at some distance from it. Never expose the patient directly to the gas flow.
- Close the cylinder valves after use, let the pressure of the pressure reducer drop by holding the flow meter open, close the flow meter and then loosen the adjusting screw of the pressure reducer (not with integrated pressure reducers).

7.3 Fittings and Application Devices

- Medical Argon supplied in its cylinder must be connected to a pressure reducer containing a country-specific gas cylinder connector and the Argon Plasma coagulator (APC) to fulfil its intended purpose.
- Before use, read the manufacturer's instructions supplied with the equipment or device associated with the Medical Argon and follow the recommendations of the physician or surgeon. In particular, check the suitability and condition of the equipment.
- The fittings and application devices used must be compatible and approved for the intended use. Use standardised fittings or hoses specifically designed for Medical Argon. Use a specific connection of type W 21,80 × 1/14 (Type 6 according to DIN 477).
- Connect the cylinder to the associated device:
 - Screw fitting: align the fittings on the regulator and the cylinder.
 - Clamp fittings: Match the positioning pins to the specific holes. Screw on the equipment by hand until it cannot be tightened any further (Note: using a spanner can damage the seal).
 - Before opening the valve, check that the connections have been fitted correctly.
 - Open the cylinder valve slowly, without forcing it or opening it fully.
 - Always close the cylinder valve after use.

- Do not disconnect the connector from the cylinder without first purging the residual compressed gas.
- Before returning for refilling, remove all accessories, hoses, etc. that were not already connected on delivery.
- Gas leak: close the valve and bleed the connection system. If the leak persists, put the cylinder outside and release the gas without attempting to patch or repair the valve. Return the defective cylinder to the manufacturer.

8. CLEANING AND DISINFECTION

- Before cleaning and disinfection, close the application device and, if necessary, release the pressure from the connected fittings.
- If exterior cleaning is required, please use only a clean cloth. The cloth can be dry or moistened with clean water. We recommend wipe disinfection of the outer surfaces.
- For disinfection wipe (not spray) cylinders with 70% solution of Isopropanol (IPA) or 70–75 % Ethanol in water or as an alternative wipe (not spray) cylinders with 0.5–1.5 % solution of H_2O_2 (hydrogen peroxide) in water. If other disinfection solutions are used, check that they are compatible with brass, plastic materials of the components (including the stickers) and the medical gas.

When the application unit is used in a hospital, the

requirements and practices of the hospital hygiene plan must be observed accordingly.

Keep alcoholic antiseptics and other flammable substances away from the electrosurgery area or other sources for ignition.

- Do not immerse the valves in water or other liquids. Do not bring any liquids into the application connections.
- Only hygienically perfect and clean cylinders (without gross contamination) may be returned.

9. MAINTENANCE, SERVICING AND DISPOSAL

→ Error – Cause – Remedy

Error	Cause	Remedy
Leakage from e.g. → Valve/cylinder connection → Filling connection → Pressure gauge on pressure reducer		Close the cylinder valve, move the containers to a ventilated area and call your Linde Gas Schweiz Service.
No gas flow, although the pressure gauge shows that the cylinder is not empty.	1. Cylinder valve is closed 2. Malfunction	1. Open the valve by turning it anticlockwise. 2. Call the Linde Gas Schweiz Service.

- No modifications may be made to the cylinder.
- Do not carry out repairs on a defective valve. Repair and maintenance may only be carried out by Linde Gas Schweiz authorised and qualified personnel.
- Unauthorised maintenance or repair will inevitably lead to the exclusion of liability.
- Linde Gas Schweiz cylinders are always to be returned to Linde Gas Schweiz and not to dispose of them yourself.

10. FURTHER INFORMATION

- Please feel free to contact Linde Gas Schweiz with any questions or anomalies.
- Argon as a medical device is classified as class IIa according to Annex VIII of the European Medical Device Regulation 2017/245.
- Argon is inert, colourless, non-corrosive, non-flammable. Argon is portion of the atmosphere (1 % by volume). The gas has the following physicochemical properties:
 - Molar mass: . . . 39,95 g/mol
 - Density, Gas 20 °C, 1 atm: 1,78 kg/m³
 - Specific Gravity, Gas (air = 1) 20 °C, 1 atm: 1,38
 - Specific Volume 20 °C, 1 atm: 0,606 m³/kg
 - 1 kg of medical gas contains at least 999,97 g Argon. The medical device contains no other ingredients.
- Further product and safety relevant information is available in the corresponding product and safety data sheet:
[Sicherheits- und Produktdatenblätter | Linde Gas Schweiz](#)
- Stay away from MR scanners. The cylinders are not compatible with magnetic fields (MRI).

→ Properties of the packaging

Article number	72030177
Cylinder	B05 Aluminium
Volume (m ³)	1,1
Pressure (bar)	200
Valve standard	W 21,80 × 1/14 (Type 6 according to DIN 477)
Valve	Residual Pressure Valve, Closure clockwise
Weight (kg)	approx. 8,0
Dimensions with cap	h × Ø (mm) approx. 610 × 140

Colour indication

Cylinder shoulder	emerald green (RAL 6001)
Cylinder body	pure white (RAL 9010)

11. LABELLING

Manufacturer

Medical Device

EC representative

Conformity mark of the medical device

0123

CE marking of the devices (incl. identification number of the notified body responsible for the conformity assessment procedures)

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Linde Gas Schweiz AG
Linde GmbH, Gases Division