**CRYOLOC®**

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**Dioxyde de carbone médical**

**Anidride carbonica medicale**

**Medical Carbon Dioxide**

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Versione 05

# ENGLISH

## Instructions for use – Information for the user

Read the entire instructions for use carefully because they contain important information. The instructions for use are an integral part of the medical device and must always be available to the user, e.g., the physician. Before working with the medical device, you must know and understand its function and handling.

The medical device CRYOLOC® consists of compressed liquefied Carbon Dioxide (CO<sub>2</sub>) gas. The liquefied gas is filled in dedicated steel or aluminium cylinders as the primary packaging, using a valve as the closure system.

Please note that a number of the points listed in this Instruction for Use may also pose a safety risk to the physician/user of CRYOLOC®. Therefore, it is particularly important that the use of CRYOLOC® must always be carried out by staff who have been trained with regard to the hazards associated with CRYOLOC®.

For detailed information, see also chapters 5, 6, and 7.

## **This instruction for use include:**

1. Intended Purpose	85
2. Intended Clinical Benefits	85
3. Intended use environment	85
4. Indication, Dosage, Method and Duration of Use	86
5. Contraindications, Warnings, Precautions, Side-Effects	88
6. Precautions and Warnings for Transport and Storage	94
7. Safe Withdrawal	95
8. Cleaning and Disinfection	100
9. Maintenance, Servicing and Disposal	102
10. Further Information	103
11. Labelling	105

## **1. INTENDED PURPOSE**

CO<sub>2</sub> is intended to be used in combination with other medical devices for the following purposes:

- As cooling/freezing media for topical cryotherapy

## **2. INTENDED CLINICAL BENEFITS**

Medical CO<sub>2</sub> is a well-established freezing/cooling medium which has a proven history to be safe and effective. The intended effect of CRYOLOC® depends on its inherent cryogenic property.

This property is used

- for freezing or ablation of cells and/or tissue
- for cooling in order to relieve pain or to decrease effusion, inflammation, or swelling
- to test the nerve response of tooth pulp suspected of being necrotic

Pharmacological effects are not intended nor expected.

## **3. INTENDED USE ENVIRONMENT**

Health care facilities.

No specific ambient conditions apply, other than those listed under section 6 «Precautions for Transport/Storage».

## **4. INDICATION, DOSAGE, METHOD AND DURATION OF USE**

### **4.1 Medical Indication**

CO<sub>2</sub> is intended to be used directly or by using probes or applicators

- as freezing media for cryosurgery
- as cooling media for topical cryotherapy like cryoanalgesia/cooling
- as a cold stimulus to test the tooth pulp health

Any indication must comply with the instructions supplied with the device connected to the gas cylinder.

### **4.2 Intended Patient Population**

Medical CO<sub>2</sub> for topical cryotherapy has no restriction regarding patient group, but for cryoanalgesia/cooling or pulp testing, patients must be able to perceive and communicate pain.

### **4.3 Dosage and Duration of Use**

For cryosurgery, cryoanalgesia/cooling, and for dental pulp testing, the duration of application, the amount of CO<sub>2</sub> used, and the flow rate depends on the therapeutic or diagnostic procedure as well as the therapeutic objective.

The physician determines the duration of application, whenever possible, according to clinical guidelines. The duration of application should be limited to the shortest possible time frame to avoid side effects.

Repeated applications depend on the intended clinical outcome, and individual patient condition and must be considered individually by the treating physicians.

#### **4.4 Intended users**

Medical personnel (e.g., physicians, dentists, nurses)

Application should only be carried out by healthcare professionals who have experience in the field of cryosurgery, cryotherapy, or endodontics and are trained in the handling of CO<sub>2</sub>.

Personnel for logistics of the cylinders (e.g., technical personnel, nurses)

No demographic or physical restrictions when working according to the Linde Gas Schweiz safety data sheet, which means that users should be trained in the handling of CO<sub>2</sub>.

#### **4.5 Method and application**

Applicators and other auxiliary devices may be required for application of Medical CO<sub>2</sub>. The type of application is determined by the applicator used to apply CO<sub>2</sub>. It should be carried out with a system that corresponds to the current state of medical technology and is compatible with (cryogenic) CO<sub>2</sub>.

Applicators are purchased by and operated under the responsibility of the users and the user should consult the user information of these devices before use.

The health care professional should adequately inform the patient about the therapy of choice, the individual cooling effect prior to the therapy application or expected diagnostic effect.

## **5. CONTRAINDICATIONS, WARNINGS, PRECAUTIONS, SIDE-EFFECTS**

- Should any incidents or side-effects occur during the use of CRYOLOC<sup>®</sup>, which are not described in chapters 5.1.3, 5.2.3, or 5.3.3, these must be reported to the manufacturer Linde Gas Schweiz.
- If you notice side-effects that are already described in chapter 5.1.3, 5.2.3, or 5.3.3, you can also report them to the manufacturer Linde Gas Schweiz.
- If a serious incident occurs in relation to CRYOLOC<sup>®</sup>, this serious incident must be reported immediately to the manufacturer and to the competent authority of the Member State in which the user and/or the patient is established.

### **5.1 Cryosurgery**

#### **5.1.1 Contraindications**

For the following indications, an alternative treatment should be chosen if possible:

- Undiagnosed skin lesions (for which pathology is required)
- Circulatory damaged areas or patients with poor circulation
- Cold/sensory sensitivity
- Previous sensitivity or adverse reaction to cryotherapy
- Dark skin or visible densely haired areas

### 5.1.2 Warnings and precautions

- Cryoablation in the airways or during upper endoscopic procedures may also dilute the inhaled gas and lead to reduction/desaturation of oxygen saturation. Patients should be monitored with pulse oximetry if Carbon Dioxide is used during tracheo/bronchoscopy or oesophagoscopy and gastroscopy, as there is a risk of dilution of the inspired gas and subsequent desaturation. 
- In the case where a larger volume of liver tissue is removed, there is a risk of cryoshock syndrome (cryosurgery of liver tissue). 
- Direct application of Carbon Dioxide, e. g., with a spray applicator, is not recommended if the skin is excessively damaged, brittle, or torn.
- Cryogenic Carbon Dioxide not only causes a freezing or cooling effect, but also leads to a space requirement due to the expansion of the gas. When Carbon Dioxide is used in body cavities, it must be ensured that the gas can flow out to avoid severe tissue damage. 
- For the patient, particular caution and eye protection should be considered when treating the facial area.
- Treatment of malignant lesions with poorly defined margins may result in recurrence.

### 5.1.3 Side-effects

When used as a freezing agent, the following side-effects are known:

- Procedure related pain and pain after treatment when warming up of the flanking tissue
- Cold pain or burning sensation
- Inflammatory reactions with swelling, oedema, redness, and pain
- Skin irritation or lesions in the treated area including frostbite, chilblains, blister formation, possibly haemorrhagic blisters
- Necrosis (in addition to the desired effect)
- Pigmentary changes and/or hair loss
- Secondary infection
- Nerve or tendon damage

### 5.1.4 Interactions

Interaction of CRYOLOC® with other medication has not been reported.

### 5.1.5 Pregnancy and lactation

There is no scientific evidence to date on the use of the product during pregnancy and breast feeding.

### 5.1.6 Roadworthiness and operation of machines

No special precautions are necessary.

## 5.2 Cryoanesthesia/Cooling

### 5.2.1 Contraindications

For the following indications, an alternative treatment should be chosen if possible:

- Circulatory damaged areas or patients with poor circulation
- Cold/sensory hypersensitivity
- Previous sensitivity or adverse reaction to cryotherapy
- Sensory impairment
- Conditions associated with blood vessel inflammation (e.g., vasculitis)
- Paroxysmal cold haemoglobinuria; cold temperatures can promote intravascular haemolysis
- Sickle cell disease; temperature changes can trigger a sickle cell crisis
- Certain autoimmune diseases may be triggered or worsened when exposed to cold temperatures. For patients suffering from autoimmune disease, immune modulating comorbidities or immune suppressive medication, the individual risk-benefit profile must be considered.

### 5.2.2 Warnings and precautions

- Avoid systemic hypothermia of the patient
- Avoid freezing! Excessive exposure may cause freezing damage (see chapter cryosurgery)
- User should take note of patient's feedback and react accordingly

- Patients must be able to perceive and communicate pain
- For the patient, particular caution and eye protection should be considered when treating the facial area
- For large-surface treatments of patients with symptomatic cardiovascular or lung disease, or blood pressure abnormalities, the individual risk-benefit profile must be considered

### 5.2.3 Side-effects

When used as a cooling agent, the following side-effect is known:

- Cold pain

### 5.2.4 Interactions

Interaction of CRYOLOC® with other medication has not been reported.

### 5.2.5 Pregnancy and lactation

There is no scientific evidence to date on the use of the product during pregnancy and breast feeding.

### 5.2.6 Roadworthiness and operation of machines

No special precautions are necessary.

## 5.3 Pulp Testing

### 5.3.1 Contraindications

For the following indications, an alternative treatment should be chosen if possible:

- Cold/sensory hypersensitivity or previous adverse reaction to cryotherapy
- Sensory impairment or nerve damage

### 5.3.2 Warnings and precautions

- Using CO<sub>2</sub> snow for cold pulp testing should not be conducted alone but along with other diagnostic test methods for an accurate diagnosis
- Patients must be able to perceive and communicate pain
- Cold tests should be applied to the tooth or teeth in question and also a control tooth until the patient reacts to the stimulus. Excessive exposure may cause freezing damage (see chapter cryosurgery)
- Caution when using a cold stimulus, on teeth that form a part of a metal bridge. The cold stimulus applied to a necrotic tooth may be felt by an adjacent tooth, as the coldness is conducted through the metal.

### 5.3.3 Side-effects

- Tissue damage of the pulp (degeneration) if tissue freezing occurs
- Fissures in the tooth enamel
- Nerve damage
- Cold pain

### 5.3.4 Interactions

Interaction of CRYOLOC® with other medication has not been reported.

### 5.3.5 Pregnancy and lactation

There is no scientific evidence to date on the use of the product during pregnancy and breast feeding.

### 5.3.6 Roadworthiness and operation of machines

No special precautions are necessary.

## 6. PRECAUTIONS AND WARNINGS FOR TRANSPORT AND STORAGE

- The carbon dioxide is supplied in pressurised cylinders. These cylinders contain both gaseous and liquid product.
- Only store and transport gas cylinders with closed valves and with the protective device provided (e.g., protective cap).
- Protect the valve from mechanical overstressing.
- During storage and transportation, restrain cylinders with more than 3.5 liters using appropriate means (chains, hooks, etc.) to keep them in a vertical position and prevent them from falling.
- Store at ambient temperature below 50°C and do not expose to sunlight or heat, burst plate of the cylinder may burst if heated.

- Store in a clean, well-ventilated room. Gas/vapours are heavier than air. They can accumulate in closed rooms, especially on the floor or in low-lying areas.
- Therefore, do not place cylinders in stairwells, corridors, passageways and recreation or consumption rooms.
- Empty and full cylinders must be labelled and stored separately from each other.  
When storing empty cylinders ensure that the valves are kept shut.
- The expiry date is printed on a batch label on the cylinder. Do not use the medical device after this date! The expiry date also applies to opened cylinders without restriction.  
The retest date according ADR requirements of the receptacle is indicated on the neck of the cylinder; the retest date of the test refers exclusively to the possibility of transport by road and does not refer to the expiry date of the product.

## **7. SAFE WITHDRAWAL**

### **7.1 General safety instructions**

- CO<sub>2</sub> may only be removed from upright cylinders. Otherwise, there is a risk of liquid carbon dioxide escaping.
- Before opening the valve, pay attention to the valve, it might be cold which can lead to cold burns.

- Opening the valve abruptly can lead to uncontrolled release of larger amounts of product and therefore introduce the risk of cryogenic burns (cold). In the event of a burn, flush with plenty of lukewarm water. Contact of solid carbon dioxide with skin or eyes leads to frostbite (cold burns) or severe eye damage.
- Using the product at a high flow rate (> 5l/min) may cause ice on the cylinders and connections.
- Ingress of liquid carbon dioxide into application device can cause serious malfunctions.
- Carbon dioxide is heavier than air. At high concentrations there is a danger of asphyxiation. Use only in rooms with sufficient ventilation to maintain CO<sub>2</sub> concentrations within the occupational exposure limits (OELs). Apply only by qualified personnel.
- In the event of gas escapes, close the leaking valve, ventilate the room well and evacuate it. Never use leaking cylinders and check that emergency measures have been taken.
- Check the cylinder before use to ensure that it is in perfect condition and suitable for the intended use.
- Remove the tamper-evident seal from the cylinder; ensure that the cylinder's fitting is clean and free of foreign bodies. Open the cylinder valve slowly until a hissing sound can be heard, then close it immediately (this procedure purges any foreign bodies). Repeat the procedure 2 to 3 times.

- In the event of anomalies, e. g., in the product application, pressure indication or labelling, do not use the product and contact the manufacturer.
- It should be noted that the gas pressure in the cylinders remains constant regardless of the level of the remaining liquid (57.3 bar at 20 °C) and thus does not allow any conclusion to be drawn about the remaining quantity. As CO<sub>2</sub> is a compressed liquefied gas, the weight will decrease quickly as soon as there is no more liquid in the cylinder. Therefore, during use, only the weight of the cylinder can give an indication of the remaining liquid level.

CO<sub>2</sub> is available in pressurised cylinders, for safe withdrawal the following must be observed:

- Before each use, the compressed gas cylinder must be checked gravimetrically for the amount of CO<sub>2</sub> present.

Safe removal of individual cylinder

The filling quantity (kg) results from the product of the filling factor (0.75) and the specific cylinder volume, stamped on the cylinder shoulder.

Example 2-litre cylinder:  $2 \text{ l} \times 0.75 \text{ kg/l} = 1.5 \text{ kg CO}_2$

Actual weight minus tare weight results in actual quantity of CO<sub>2</sub> available.

Safe removal of bundles

Before use, check the bundle for residual pressure.

Bundle à 12 single cylinders. Quantity approx. 450 kg.

Ensure that always a reserve cylinder / bundle is available.

- Do not lift cylinders by their valves.
- Never squeeze a cylinder into a holder, into which it is difficult to fit.
- Misuse of the pressure containers, filling by the user or by third parties and decanting into other cylinders are not permitted. Do not transfer compressed gas from one cylinder to another.
- The cylinders must not be completely emptied; they must be returned with a residual pressure of at least 3 bar. This residual pressure is necessary to protect the gas cylinder from contamination and liquid ingress or to prevent corrosion due to moisture.

## **7.2 Operating CO<sub>2</sub> cylinders with dip tube**

- Gas cylinders with a dip tube are marked with a label outlined in red. They may only be used when liquid carbon dioxide is to be withdrawn.
- Do not use cylinders with dip tubes for insufflation, due to risk of cold burns related to leakage of liquefied product)!
- Never open the CO<sub>2</sub> cylinder equipped with a dip tube without using a connector (risk of cold burns from escaping liquefied product).
- Do not use a pressure reducer for cylinders with a dip tube, as the liquid carbon dioxide cools down considerably due to expansion and a mixture of CO<sub>2</sub> gas and CO<sub>2</sub> snow is formed and could clog the pressure reducer.

- Carbon dioxide may only be taken from upright cylinders, so that the opening of the dip tube remains below the CO<sub>2</sub> liquid level.
- The valve of the cylinder must be opened slowly (turn anticlockwise), otherwise there is a risk of cold burn!
- The withdrawal is carried out with full cylinder pressure. The application device must be appropriately pressure-resistant and resistant to liquid carbon dioxide.

### **7.3 Fittings and application devices**

- Please verify that your cylinder is intended to provide liquid CO<sub>2</sub>.
- CO<sub>2</sub> supplied in its cylinder must be connected to one or more other medical devices (usually a device for cryosurgery, a device for cryotherapy or a device for Cold Pulp Testing (Odontotest)) to fulfil its intended purpose. The choice of suitable device depends on the application and is the responsibility of the user.
- Before use, read the manufacturer's instructions supplied with the equipment or device associated with the medical CO<sub>2</sub> and follow the recommendations of the physician or surgeon. In particular, check the suitability and condition of the equipment.
- The fittings and application devices used must be compatible and approved for the intended use. Use standardised fittings or hoses specifically designed for medical CO<sub>2</sub>. Use a specific connection of type DIN 477 No. 6 (W21,8×1/14")

- Connect the cylinder to the associated device:
  - Screw fitting: align the fittings on the regulator and the cylinder.
  - Clamp fittings: Match the positioning pins to the specific holes.
  - Screw on the equipment by hand until it cannot be tightened any further (in some cases, using a spanner can damage the seal).
  - Before opening the valve, check that the connections have been fitted correctly.
  - Open the cylinder valve slowly, without forcing it or opening it fully.
  - Always close the cylinder valve after use
  - Do not disconnect the connector from the cylinder without first purging the residual compressed gas.
- Before returning for refilling, remove all accessories, hoses, etc. that were not already connected on delivery.
- Gas leak: close the valve and bleed the connection system. If the leak persists, put the cylinder outside and release the gas without attempting to patch or repair the valve. Return the defective cylinder to the manufacturer.

## **8. CLEANING AND DISINFECTION**

- Before cleaning and disinfection, close the application device and, if necessary, release the pressure from the connected fittings.

- If exterior cleaning is required, please use only clean cloth. The cloth can be dry or moistened with clean water. We recommend wiping disinfection of the outer surfaces.
- For disinfection wipe (not spray) cylinders with 70 % solution of Isopropanol (IPA) or 70–75 % Ethanol in water or as an alternative wipe (not spray) cylinders with 0.5–1.5 % solution of H<sub>2</sub>O<sub>2</sub> (hydrogen peroxide) in water. If other disinfection solutions are used, check that they are compatible with brass, plastic materials of the components (including the stickers) and the medical gas.  
When the application unit is used in a hospital, the requirements and practices of the hospital hygiene plan must be observed accordingly.
- Do not immerse the valves in water or other liquids. Do not bring any liquids into the application connections.
- Only hygienically perfect and clean cylinders (without gross contamination) may be returned.

## 9. MAINTENANCE, SERVICING AND DISPOSAL

→ Error - Cause - Remedy

Error	Cause	Remedy
Leakage from e. g. → Valve/cylinder connection → Filling connection → Pressure gauge on pressure reducer → Bundle screw connection		Close the cylinder valve/ bundle valve, move the containers to a ventilated area and call your Linde Gas Schweiz Service.
No gas flow, although the pressure gauge shows that the cylinder is not empty.	1. Cylinder valve is closed. 2. Malfunction	1. Open the valve by turning it anticlockwise. 2. Call the Linde Gas Schweiz Service.

- No modifications may be made to the cylinder.
- Do not carry out repairs on a defective valve. Repair and maintenance may only be carried out by Linde Gas Schweiz authorised and qualified personnel.
- Unauthorised maintenance or repair will inevitably lead to the exclusion of liability.
- Linde Gas Schweiz cylinders are always to be returned to Linde Gas Schweiz and not to dispose of them yourself.

## 10. FURTHER INFORMATION

- Please feel free to contact Linde Gas Schweiz with any questions or anomalies.
- Carbon Dioxide as a medical device is classified as class IIa according to Annex VIII of the European Medical Device Regulation (EU) 2017/745.
- Carbon dioxide is a very stable, non-flammable and inert, colourless and odourless gas that is heavier than air and has the following physicochemical properties:
  - Molar mass: 44,010 g/mol
  - Sublimation point at 1.013 bar: 194.65 K (-78,5 °C)
  - Vapour pressure at 20 °C: 57.3 bar
  - 1 kg of medical gas contains as active ingredient component: carbon dioxide, at least 995 g.  
The medical device contains no other ingredients.
- Further product and safety relevant information is available in the corresponding product and safety data sheet:  
[Sicherheits- und Produktdatenblätter | Linde Gas Schweiz](#)
- Stay away from MR scanners. The cylinders are not compatible with magnetic fields (MRI).
- The cylinder is equipped with a valve, with rupture disc.

→ Properties of the packaging

Material number	3750172VT2	3750173VT2	3750112VT2	37501811VT2	37505883VT2
Cylinder	Steel	Steel	Steel	Steel/ Aluminium	Steel
Cylinder Volume (Liter)	2	3,5	10	40	600
Valve standard	DIN 477 No. 6 (W21,8 x 1/14")	DIN 477 No. 6 (W21,8 x 1/14")	DIN 477 No. 6 (W21,8 x 1/14")	DIN 477 No. 6 (W21,8 x 1/14")	DIN 477 No. 6 (W21,8 x 1/14")
Valve	RPN clockwise	RPN clockwise	RPN clockwise	RPN clockwise	RPN clockwise
Dip tube	Yes	Yes	Yes	Yes	Yes
Weight (kg)	5	7,5	26	76	1500
Dimensions with cap h x Ø (mm)	415 x 102	650 x 100	930 x 140	1610 x 205	Bundle Height 196 cm Length 100 cm Width t 75 cm

RPN – Residual Pressure Valve      Colour indication:      Cylinder shoulder dusty grey (RAL 7037)  
Cylinder body pure white (RAL 9010)

## 11. LABELLING

Manufacturer

Medical Device

EU representative

Conformity mark of the medical device

CE marking of the devices (incl. identification number of the notified body responsible for the conformity assessment procedures)

Catalogue number, Article number

Unique Device Identifier

Batch number

Use-by date, Expiry date

Caution, warning notice

Electronic instruction for use  
Read the electronic instructions for use before  
use at

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For printed instructions for use,  
call 0844 800 300.

Material Number	Material Name	UDI-DI
3750172	LAPAROX® CO <sub>2</sub> MED 2 L	111756893460
3750173	LAPAROX® CO <sub>2</sub> MED 3,5 L	111756894026
3750173V09	LAPAROX® CO <sub>2</sub> MED 3,5 L PI	111756895716
3750112	LAPAROX® CO <sub>2</sub> MED 10 L	111756896379
37501811	LAPAROX® CO <sub>2</sub> MED 40 L	111756898635
37505883	LAPAROX® CO <sub>2</sub> MED 600 L	111756899201

Linde Gas Schweiz AG  
Industriepark 10  
CH-6252 Dagmersellen  
Telefon 0844 800 300, Fax 0844 800 301  
contact.lg.ch@linde.com, www.linde.ch

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Linde GmbH, Gases Division