



GLOOR JETFLOW (lynkobling).



GLOOR JETFLOW (gevindtilslutning).

Produktbeskrivelse Flowvælger GLOOR JETFLOW findes som standard til medicinsk oxygen. Flowvælgerne leveres med lynkobling som passer i gasudtag for medicinsk oxygen og med gevindtilslutning til montering i ambulance.

Vedligeholdelse Enheden skal gennemgå fuld service hvert sjette år. Service må kun udføres af autoriseret personale.

Tekniske data	Klassificering	Enhedsklasse for "EU guideline" 93/42 EEC Annex IX: Klasse IIb
	Standard	EN ISO 15002
	Konstruktionstype	Flowkontrol med flowvælger
	Forsyningsstryk	4,5 + 0,2 bar
	Nøjagtighed	+/- 20 %. Kalibreret for 450 kPa (4,5 bar) Flowvælgerne er modtrykstilkomponerede

Flowomåder	GLOOR JETFLOW 0-5 l/min	0	0,5	1,0	1,5	2,0	2,5	3,0	4,0	5,0
	GLOOR JETFLOW 0-15 l/min	0	1,0	2,0	3,0	4,0	6,0	8,0	12,0	15,0
	GLOOR JETFLOW 0-31 l/min	0	3,0	5,0	9,0	12,0	15,0	20,0	26,0	31,00
	GLOOR JETFLOW 1-31 l/min	0	1,0	3,0	5,0	10,0	15,0	20,0	26,0	31,0

Artikelnummer	Art. nr.	Betegnelse
	320 781	GLOOR JETFLOW 0-5 l/min med lynkobling
	320 768	GLOOR JETFLOW 0-15 l/min med lynkobling
	320 769	GLOOR JETFLOW 0-31 l/min med lynkobling
	331 389	GLOOR JETFLOW 1-31 l/min med lynkobling
	320 780	GLOOR JETFLOW 0-15 l/min til montering i ambulance
	320 783	GLOOR JETFLOW 0-31 l/min til montering i ambulance